

COMPLIMENTARY WEBINAR

Know the Pneumonia: Urinary Antigen Testing in Patient Care

Thursday, November 9, 2023

1:00 PM – 2:00 PM ET



JANET E. STOUT, PHD

Executive V.P. and Founder
Special Pathogens Laboratory
Pittsburgh, PA



NORMAN MOORE, PHD

Director, Medical Affairs
Infectious Diseases, North America
Abbott
Scarborough, ME

The speakers are presenting on behalf of Abbott.

The information presented is consistent with applicable FDA guidelines.

This program does not provide continuing medical education (CME) credits.

© 2023 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. COL-22301-01 10/23

Sponsored by:



Joining the Teleconference


If you closed the 'Phone' call in box when you joined:



← Click 'Webinar Info' for Dial-in information
(top left of your screen)

Webinar communication tools

Show **Chat** panel



Other panels:

Q&A

Notes

Webinar communication tools

Chat →

Chat [close]

Type here →

To: [dropdown]
Enter chat message here [send]

Webinar communication tools

Chat →

Type here →

Chat [Close] [Share]

To: [Dropdown]

Enter chat message here [Send]

Q&A →

Type here →

Q & A [Close] [Share]

All (0)

Select a question and then type your answer here. There's a 512-character limit.

[Send] [Send Privately]

Webinar communication tools

Chat →

Chat [Close] [X]

To: [Dropdown]

Enter chat message here [Send]

Type here →

Q&A →

Q & A [Close] [X]

All (0)

Select a question and then type your answer here. There's a 512-character limit.

[Send] [Send Privately]

Type here →

CC [Close] ← Closed Captioning

Available CE Credit

P.A.C.E.[®]

Florida laboratory CE

Certificate of Attendance

To Receive Certificate of Attendance

After today's webinar:

- A certificate of attendance available for all attendees
- Evaluation form will appear automatically
- Must complete Eval to receive Certificate link via email
- For groups: Those logged in will receive Email from messenger@webex.com with link to evaluation.
Forward email to colleagues who attended with you!!!
- Double-check email address

Joined Using a Mobile Device?

- Evaluation won't appear automatically, but...
- Watch for email with link to evaluation!



Recording

Within a few days following today's event, visit
<https://www.whitehatcom.com/abbott>



Know the Pneumonia: Urinary Antigen Testing in Patient Care

Live Event: Thursday, November 9, 2023 | 1:00 - 2:00 PM Eastern Time

P.A.C.E.® credit available until November 9, 2024

Florida Laboratory CE Credit available

RECORDING

SLIDES

Join this session for insights on improving the quality of care of patients with pneumonia. Experts will discuss the impact of this deadly infection and current gaps in identifying *Legionella* and *Streptococcus (S.) pneumoniae*. Practical applications of integrating urinary antigen testing (UAT) in patient scenarios will be explored, including an evidence-based review of mortality and antibiotic stewardship improvements achieved with UAT.

Presenter:



Janet E. Stout, PhD

**Executive V.P. and Founder
Special Pathogens Laboratory
Pittsburgh, PA**

Moderator/Speaker:



Norman Moore, PhD

**Director, Infectious Diseases,
Medical Affairs, Rapid Diagnostics
Abbott
Scarborough, ME**

The webinar will:

- Evaluate the impact of pneumonia and risks for severe illness
- Examine gaps and potential delays in the diagnosis of pneumonia and quality measures to enhance patient care
- Review the identification of *Legionella* and *S. pneumoniae* with UAT
- Discuss practical examples where UAT testing supported improved patient outcomes through timely identification



MODERATOR AND SPEAKER

NORMAN MOORE, PHD

Director, Infectious Diseases
Scientific Affairs
Rapid Diagnostics, Abbott

Disclosures

Employed by and speaking on behalf of Abbott

Objectives


- Evaluate the impact of pneumonia and risks for severe illness
- Examine gaps and potential delays in the diagnosis of pneumonia and quality measures to enhance patient care
- Review the identification of *Legionella* and *S. pneumoniae* with UAT
- Discuss practical examples where UAT testing supported improved patient outcomes through timely identification

Background on Urinary Antigen Testing

The background of the slide features a gradient from dark purple on the left to bright blue on the right. Overlaid on this gradient are several blurred, light blue and white shapes that resemble microscopic organisms, such as bacteria or cells, scattered across the right side of the frame.

What was the need for Urinary Antigen Testing?




- High mortality
- No practical or timely detection method
 - Blood culture
 - X-ray
- Poor sputum specimen



- Urinary Antigen Test technology developed to diagnose *Legionella*
- Urine provided ease of collection and faster detection
- Associated with reduction in *Legionella* mortality¹
- Catalyst for development of UAT for *S. pneumoniae* detection

1. Benin AL, et al. Trends in Legionnaires Disease, 1980–1998: Declining Mortality and New Patterns of Diagnosis. *Clinical Infect Dis*. 2002; 35:1039–46.

Examples of Urinary Antigen Testing





JANET E. STOUT, PHD

Executive V.P. and Founder
Special Pathogens Laboratory
Pittsburgh, PA

Update on Legionellosis: A Not So Atypical Pneumonia

Janet E. Stout, PhD

President, Special Pathogens Technology

Founder, Special Pathogens Laboratory

Disclaimers

- Abbott is the program sponsor, the content of the presentation is consistent with applicable governing regulatory body requirements, and I was chosen by Abbott and am presenting the program material on Abbott's behalf.
- No conflicts of interest to disclose

Infectious Disease Microbiologist



- Microbiologist with more than 30 years studying Legionnaires' disease

Infectious Disease Microbiologist



- Microbiologist with more than 30 years studying Legionnaires' disease
 - I'm a Legionellologist

Infectious Disease Microbiologist



- Microbiologist with more than 30 years studying Legionnaires' disease
 - I'm a Legionellologist
- Mission: What I've learned can help you address risks from *Legionella*

Objectives

Objectives

- What is Legionnaires' disease? It's pneumonia and more

Objectives

- What is Legionnaires' disease? It's pneumonia and more
- Why is Legionnaires' disease under-reported?

Objectives

- What is Legionnaires' disease? It's pneumonia and more
- Why is Legionnaires' disease under-reported?
- Making the diagnosis of Legionnaires' disease

Objectives

- What is Legionnaires' disease? It's pneumonia and more
- Why is Legionnaires' disease under-reported?
- Making the diagnosis of Legionnaires' disease
- Knowing what's in your water can inform diagnostic testing

Objectives

- What is Legionnaires' disease? It's pneumonia and more
- Why is Legionnaires' disease under-reported?
- Making the diagnosis of Legionnaires' disease
- Knowing what's in your water can inform diagnostic testing
- Some tricks of the trade

What Is Legionellosis?

Legionnaires' disease (LD) is a multi-system illness, with pneumonia, caused by *Legionella* species

What Is Legionellosis?

Legionnaires' disease (LD)

is a multi-system illness, with pneumonia, caused by *Legionella* species

Pontiac fever

is a self-limited flu-like illness, without pneumonia, that is associated with *Legionella* species ... Resolves without treatment

Council of State and Territorial Epidemiologists (CSTE)

Adopted in June 2019 and put into effect January 1, 2020



Addition - Extrapulmonary Legionellosis

- Disease at sites outside of lung (endocarditis, wound infection, joint infection)
- Now three clinically and epidemiologically distinct forms of Legionellosis.

What Are Epidemiologic Risk Factors for Legionnaires' Disease?

- Recent travel with an overnight stay outside of the home
- Exposure to whirlpool spas and hot tubs
- Recent repairs or maintenance work on domestic plumbing
- Renal or hepatic failure
- Diabetes
- Systemic malignancy
- Smoking
- Immune system disorders
- Aged more than 50 years old

Resources

CDC

Sources: US Centers for Disease Control and Prevention. *Legionella* (Legionnaires' Disease and Pontiac Fever): Disease Specifics. (Reviewed: Mar 25, 2021.) Accessed Jul 23, 2023. <https://www.cdc.gov/legionella/clinicians/disease-specifics.html>

US Centers for Disease Control and Prevention. *Legionella* (Legionnaires' Disease and Pontiac Fever): Diagnosis, Treatment & Prevention. (Reviewed: Mar 25, 2021.) Accessed Jul 23, 2023. <https://www.cdc.gov/legionella/clinicians/diagnostic-testing.html>

Viasus D, Gaia V, Manzur-Barbur C, Carratalà J. Legionnaires' Disease: Update on Diagnosis and Treatment. *Infect Dis Ther.* 2022 Jun;11(3):973-986. doi: 10.1007/s40121-022-00635-7. Epub 2022 May 3. PMID: 35505000; PMCID: PMC9124264.

Clinical Notes

Clinical Notes

Elderly – signs not typical

Clinical Notes

Elderly – signs not typical

40%–50% of Legionnaires' disease patients
have neurological symptoms

Clinical Notes

Elderly – signs not typical

40%–50% of Legionnaires' disease patients
have neurological symptoms

Some with encephalopathy (brain inflammation) triggered by
antibodies that target brain cells (www.legionella.org)

Clinical Notes

Elderly – signs not typical

40%–50% of Legionnaires' disease patients
have neurological symptoms

Some with encephalopathy (brain inflammation) triggered by
antibodies that target brain cells (www.legionella.org)

Anti-tumor necrosis factor (TNF)-alpha treatment

Legionnaires' Disease

In the U.S. approximately 1 million adults are diagnosed with community-acquired pneumonia requiring hospitalization annually

Legionnaires' Disease

In the U.S. approximately 1 million adults are diagnosed with community-acquired pneumonia requiring hospitalization annually

2 – 5% are caused by *Legionella*,
at least 30,000 cases/year (minimum)

Legionella and Healthcare Costs

Estimate of Burden and Direct Healthcare Cost of Infectious Waterborne Disease in the United States

Sarah A. Collier, Li Deng, Elizabeth A. Adam, Katharine M. Benedict, Elizabeth M. Beshearse,
Anna J. Blackstock, Beau B. Bruce, Gordana Derado, Chris Edens, Kathleen E. Fullerton,
Julia W. Gargano, Aimee L. Geissler, Aron J. Hall, Arie H. Havelaar, Vincent R. Hill, Robert M. Hoekstra,
Sujan C. Reddy, Elaine Scallan, Erin K. Stokes, Jonathan S. Yoder, Michael J. Beach

Emerging Infectious Diseases • www.cdc.gov/eid • Vol. 27, No. 1, January 2021

Results

- 7.15 million waterborne illnesses occur annually
- 6,630 deaths
- \$3.33 billion in direct healthcare costs
- \$2.39 billion (72%) of those costs due to hospitalizations and deaths from:
 - Nontuberculous mycobacteria, *Pseudomonas* and *Legionella*

Legionella and Disease Burden

Estimating Waterborne Infectious Disease Burden by Exposure Route, United States, 2014

Megan E. Gerdes,¹ Shanna Miko,¹ Jasen M. Kunz, Elizabeth J. Hannapel, Michele C. Hlavsa, Michael J. Hughes, Matthew J. Stuckey, Louise K. Francois Watkins, Jennifer R. Cope, Jonathan S. Yoder, Vincent R. Hill, Sarah A. Collier

Emerging Infectious Diseases • www.cdc.gov/eid • Vol. 29, No. 7, July 2023

Legionella A Major Contribution

- Waterborne infection types accounted for 72% of \$3.33 billion total costs
 - Legionnaires' disease for 12%
 - Pseudomonas pneumonia for 14%
 - NTM infection was responsible for 46%

Mortality and Healthcare-Acquired Legionnaires' Disease



“Mortality for Legionnaires' disease continues to be high – as high as 25 percent for cases acquired in the healthcare facility.”

Pierre DM, Baron J, Yu VL, Stout JE. Diagnostic testing for Legionnaires' disease. *Ann Clin Microbiol Antimicrob.* 2017 Aug 29;16(1):59. doi: 10.1186/s12941-017-0229-6. PMID: 28851372; PMCID: PMC5576257.

Diagnosis of Legionnaires' Disease

- Adopt a proactive approach to diagnosis of pneumonias and anticipate the possibility of Legionnaires' disease.

Pierre et al. *Ann Clin Microbiol Antimicrob* (2017) 16:59
DOI 10.1186/s12941-017-0229-6

Annals of Clinical Microbiology
and Antimicrobials

REVIEW

Open Access



Diagnostic testing for Legionnaires' disease

David M. Pierre¹, Julianne Baron^{1,2}, Victor L. Yu^{1,3*} and Janet E. Stout^{1,2}

Something To Remember

Early Diagnosis = Better Outcome

Delayed diagnosis = Increased Mortality

Delayed Diagnosis = Delayed Therapy

286

Article

Vol. 15, No. 4

Eur. J. Clin. Microbiol. Infect. Dis., 1996, 15: 286–290

Delay in Appropriate Therapy of *Legionella* Pneumonia Associated with Increased Mortality

C.H. Heath^{1,2*}, D.I. Grove¹, D.F.M. Looke^{1,3}

Early Diagnosis With UAT = Better Outcome

RESEARCH

Legionnaires' Disease Outbreak in Murcia, Spain

**Ana García-Fulgueiras,* Carmen Navarro,* Daniel Fenoll,† José García,* Paulino González-Diego,*
Teresa Jiménez-Buñuales,* Miguel Rodríguez,* Rosa Lopez,* Francisco Pacheco,* Joaquín Ruiz,‡
Manuel Segovia,§ Beatriz Baladrón,¶ and Carmen Pelaz¶**

Emerging Infectious Diseases • Vol. 9, No. 8, August 2003 915

Murcia, Spain Very Low Case-Fatality

- 449 confirmed cases and an estimated total number of cases of 650.
- The reported case-fatality rate (1%) is much lower than those observed in other community outbreaks

Murcia, Spain Very Low Case-Fatality

- 449 confirmed cases and an estimated total number of cases of 650.
- The reported case-fatality rate (1%) is much lower than those observed in other community outbreaks
- This rate can be attributed, at least partially, to the quick detection of the outbreak, early diagnosis of the disease with urine antigen testing, and appropriate treatment of patients.

Something To Remember

Legionella is the problem
you don't think you have until you have it

Be Proactive!

Diagnosis: Many Cases Missed

Diagnostic tests for *Legionella* not routine –
often not done

Diagnosis: Many Cases Missed

Diagnostic tests for *Legionella* not routine – often not done

Many studies have demonstrated under reporting/missed diagnosis

Under Reporting Due to Missed Diagnosis

Hollenbeck et al. *BMC Infectious Diseases* 2011, **11**:237
<http://www.biomedcentral.com/1471-2334/11/237>



RESEARCH ARTICLE

Open Access

How often is a work-up for *Legionella* pursued in patients with pneumonia? A retrospective study

Brian Hollenbeck¹, Irene Dupont² and Leonard A Mermel^{2,3*}

Abstract

Background: It is unclear how often patients with pneumonia are assessed for *Legionella* in endemic areas. Additionally, the sensitivity of the IDSA/ATS criteria for recommended *Legionella* testing is undefined.

Methods: We performed a single-center, retrospective study of patients diagnosed with *Legionella* pneumonia at our hospital to determine: 1) how often *Legionella* diagnostic testing is obtained on patients with pneumonia at the time of hospitalization or when pneumonia developed during hospitalization; and 2) how often patient's with *Legionella* pneumonia met at least one of the five criteria in the IDSA/ATS guidelines recommending a work-up for

Results for 37 Cases

41% of *Legionella* cases were missed when following current IDSA-ATS recommendations for *Legionella* testing

IDSA – Infectious Disease Society of America
ATS – American Thoracic Society

Hollenbeck, B., Dupont, I., & Mermel, L. A. (2011).

Clinical Cultures

Providers should be educated to maintain a high index of suspicion for diagnosis of healthcare-associated Legionellosis

Clinical Cultures

Providers should be educated to maintain a high index of suspicion for diagnosis of healthcare-associated Legionellosis

Should perform cultures for *Legionella* of appropriate respiratory specimen **and** urinary antigen test

Legionella Urinary Antigen Tests (UAT)

Enzyme Immunoassay (EIA)

Immunochromatographic Test (ICT)



<https://www.globalpointofcare.abbott/us/en/product-details/binaxnow-legionella-urinary-antigen-eia.html>

<https://www.globalpointofcare.abbott/us/en/product-details/binaxnow-legionella-us.html>

<https://www.meridianbioscience.com/diagnostics/disease-areas/respiratory/legionella/tru-legionella/>

<https://immuvue.com>

Diagnostic Methods: UAT Rules!

Urine antigen tests (UAT) confirmed
97% of U.S. cases

Diagnostic Methods: UAT Rules!

Urine antigen tests (UAT) confirmed
97% of U.S. cases

Less than 10% of cases confirmed by culture –
mostly because culture was not ordered

Physicians Should Order Both Culture and Urine Antigen!

What Clinicians Need to Know about

LEGIONNAIRES' DISEASE

Legionnaires' disease is a sometimes fatal form of pneumonia that is on the rise in the United States. Unfortunately, this disease is also underrecognized and underdiagnosed. Clinicians are in a unique position to make sure cases are detected, allowing rapid investigation by public health officials and prevention of additional cases.

Diagnosis and Testing

Clinical features of Legionnaires' disease include cough, fever, and radiographic pneumonia. Signs and symptoms for Legionnaires' disease are similar to pneumonia caused by other pathogens; the only way to tell if a pneumonia patient has Legionnaires' disease is by getting a specific diagnostic test. Indications that warrant testing include:

- Patients who have failed outpatient antibiotic therapy for community-acquired pneumonia
- Patients with severe pneumonia, in particular those requiring intensive care
- Immunocompromised patients with pneumonia*
- Patients with a travel history (patients who have traveled away from their home within 10 days before the onset of illness)

Order both a culture of a lower respiratory specimen and a urinary antigen test when testing patients for *Legionella*.

Legionnaires' Disease Is On the Rise
2000–2018*



Physicians Should Order Both Culture and Urine Antigen!

What Clinicians Need to Know about

LEGIONNAIRES' DISEASE

Legionnaires' disease is a sometimes fatal form of pneumonia that is on the rise in the United States. Unfortunately, this disease is also underrecognized and underdiagnosed. Clinicians are in a unique position to make sure cases are detected, allowing rapid investigation by public health officials and prevention of additional cases.

Diagnosis and Testing

Clinical features of Legionnaires' disease include cough, fever, and radiographic pneumonia. Signs and symptoms for Legionnaires' disease are similar to pneumonia caused by other pathogens; the only way to tell if a pneumonia patient has Legionnaires' disease is by getting a specific diagnostic test. Indications that warrant testing include:

- Patients who have failed outpatient antibiotic therapy for community-acquired pneumonia
- Patients with severe pneumonia, in particular those requiring intensive care
- Immunocompromised patients with pneumonia*
- Patients with a travel history (patients who have traveled away from their home within 10 days before the onset of illness)

Order both a culture of a lower respiratory specimen and a urinary antigen test when testing patients for *Legionella*.

Legionnaires' Disease Is On the Rise
2000–2018*



Sensitivity and Specificity of Diagnostic Tests

Sensitivity varies depending on the quality and timing of clinical specimen collection, as well as technical skill of the laboratory worker performing the test. The table below provides general ranges for the sensitivity and specificity of each diagnostic test.

Test	Sensitivity (%)	Specificity (%)
Culture	20–80	100
Urinary antigen for <i>L. pneumophila</i> serogroup ¹ (Lp1)	70–100	95–100
Polymerase Chain Reaction (PCR) ²	95–99	>99
Direct Fluorescent Antibody (DFA) Stain	25–75	>95
Paired serology ³	80–90	>99



¹ Cross reactions with other species and serogroups have been documented.

²Avni T, Bieber A, Green H, et al. [Diagnostic accuracy of PCR alone and compared to urinary antigen testing for detection of *Legionella* spp.: A](https://www.cdc.gov/legionella/clinicians/diagnostic-testing.html)

Urine Antigen Sensitivity/Specificity

Duration

Excretion can be long
(weeks and sometimes months
in immuno-compromised)

Urine Antigen Sensitivity/Specificity

Duration

Excretion can be long
(weeks and sometimes months
in immuno-compromised)

False
negative
results

In confirmed cases

Urine Antigen Sensitivity/Specificity

Duration

Excretion can be long
(weeks and sometimes months
in immuno-compromised)

False
negative
results

In confirmed cases


False
positive
results

Non-specific binding in EIA format

Proteins in urine

Early Detection And Duration of Urine Antigen Positivity

Sequential Urines from 21 Culture Positive Patients



*n = no. patients

Stout - Personal communication

Why Combine Urine Antigen With Culture?


Urine antigen specific for
Legionella pneumophila (Lp), serogroup 1 only

Why Combine Urine Antigen With Culture?

Urine antigen specific for
Legionella pneumophila (Lp), serogroup 1 only

Legionnaires' disease has been caused by
other serogroups and species of *Legionella*

Legionella pneumophila, serogroup 1: Most virulent and cause of most infections



Know What's In Your Water

The urine antigen test can help detect cases of healthcare-acquired Legionnaires' disease

If used to screen for healthcare-acquired LD, you'd better know what's in your water!

- If Lp-6 in the water, diagnosis will be missed

Tips

- Create an electronic ward order to prompt physicians to order both *Legionella* urine antigen and culture.

Tips

- Create an electronic ward order to prompt physicians to order both *Legionella* urine antigen and culture.
- If you suspect that a bloody urine or proteins in the urine have caused a false positive EIA urine antigen test, follow with an ICT test.

Tips

- Create an electronic ward order to prompt physicians to order both *Legionella* urine antigen and culture.
- If you suspect that a bloody urine or proteins in the urine have caused a false positive EIA urine antigen test, follow with an ICT test.
- Do not discard sputum if epithelial cells present


Tips

- Create an electronic ward order to prompt physicians to order both *Legionella* urine antigen and culture.
- If you suspect that a bloody urine or proteins in the urine have caused a false positive EIA urine antigen test, follow with an ICT test.
- Do not discard sputum if epithelial cells present.
- If urine antigen test is positive, retrieve the sputum that was ordered for routine work-up.

Save the Sputum!



End Legionnaires' Disease





Know what's in your water so you can respond, not react



THANK YOU!

Janet E. Stout, PhD
President, Special Pathogens Technology
Founder, Special Pathogens Laboratory

E: jstout@specialpathtech.com

W: SpecialPathogensTechnology.com



NORMAN MOORE, PHD

Director, Infectious Diseases
Scientific Affairs
Rapid Diagnostics, Abbott


PNEUMONIA

S. pneumoniae



S. pneumoniae

Most common pathogen for community-acquired pneumonia (CAP)¹




Most common secondary infection with influenza diagnosis³


1. Dion CF, Ashurst JV. Streptococcus Pneumoniae. [Updated 2022 Aug 8]. In: StatPearls [Internet]. Treasure Island (FL): StatPearls Publishing; 2022 Jan.
2. CDC. Pneumococcal Disease. Fast Facts. <https://www.cdc.gov/pneumococcal/about/facts.html>, updated January 27, 2022.
3. Morris et al. Secondary Bacterial Infections Associated with Influenza Pandemics. *Frontiers in Microbiology*. 2017; 8:1041.

Etiology of Community-Acquired Pneumonia


A Specific Pathogens Detected




Etiology of Community-Acquired Pneumonia



Streptococcus pneumoniae Drug Resistance



Drug-resistant *Streptococcus pneumoniae*



*Unable to compare data with 2019 report estimates, see [Methods](#) for details.

Symptoms of Pneumococcal Pneumonia Can Vary

Pneumococcal disease can include many different [types of infections](#). Symptoms depend on the part of the body that is infected. Most pneumococcal infections are mild. However, some can be deadly or result in long-term problems.

Symptoms



Headache



Stiff neck



Fever
or chills



Confusion



Difficulty
breathing



Sensitivity
to lights



Ear pain



Cough



Chest pain



Sore throat

“Older adults with pneumococcal pneumonia may experience confusion or low alertness, rather than the more common symptoms listed above.” CDC

Diagnostic Tests for *S. pneumoniae*

METHODOLOGY	TURNAROUND TIME	SENSITIVITY	SPECIFICITY	SAMPLE TYPE	COMPONENT DETECTED
UAT¹	15 minutes	86%*	94%*	urine	antigen
Sputum Gram Stain⁶	15 minutes	15% - 100%	11% - 100%	sputum	organism
Blood Culture	24 - 48 hours	10% - 30% ²	N/A	blood	organism
Sputum Culture	24 - 48 hours	29% - 94% ³	66% ⁴ - 94% ⁵	sputum	organism

* Sensitivity and specificity data are retrospective for urine only.

1. BinaxNOW™ *S. pneumoniae* Urinary Antigen Card Package Insert.
2. Schrag SJ, et al. Resistant Pneumococcal Infections, WHO/CDS/CSR/DRS/2001.6.
3. Musher D, et al. Diagnostic Value of Microscopic Examination of Gram-Stained Sputum and Sputum Cultures Inpatients with Bacteremic Pneumococcal Pneumonia; CID: 2004:39.
4. Stralin K, et al. Etiologic Diagnosis of Adult Bacterial Pneumonia by Culture and PCR Applied to Respiratory Tract Samples, J Clin Micro, Feb. 2006, 643-645.
5. Garcia-Vazquez E, et al. Assessment of the Usefulness of Sputum Culture for Diagnosis of Community-Acquired Pneumonia Using the PORT Predictive Scoring System, Arch Inter Med/Vol. 164, Sept. 13, 2004, 1807-1811.
6. Reed, W, et al. Sputum Gram's Stain in Community Acquired Pneumococcal Pneumonia – A Meta-analysis; West J. Med 1996; 165:197-204.

ATS/IDSA CAP Clinical Practice Guidelines, 2019

In adults with CAP
***Legionella* and pneumococcal urinary antigen testing should be performed at the time of diagnosis**

C Conditional

Recommendation

***S. pneumoniae* and *Legionella* antigen testing**

- **Adults with severe CAP**
- Where indicated by **epidemiological factors (*Legionella*)**
 - i.e., known outbreaks or recent travel

Minor criteria

Respiratory rate ≥ 30 breaths/min
 $Pa_{O_2}/F_{I_{O_2}}$ ratio ≤ 250
Multilobar infiltrates
Confusion/disorientation
Uremia (blood urea nitrogen level ≥ 20 mg/dl)
Leukopenia* (white blood cell count $< 4,000$ cells/ μ l)
Thrombocytopenia (platelet count $< 100,000$ / μ l)
Hypothermia (core temperature $< 36^\circ\text{C}$)
Hypotension requiring aggressive fluid resuscitation

Major criteria

Septic shock with need for vasopressors
Respiratory failure requiring mechanical ventilation

Importance of Flu Testing During Respiratory Season

S. pneumoniae: A Secondary Complication to Influenza

- 2009 pandemic Influenza A (H1N1) & Spanish flu 1918
 - Many deaths were attributed to the flu combined with the secondary complication of pneumonia¹
- Testing for *S. pneumoniae* and Influenza helps inform the appropriate use of antibiotic therapy²
 - Is it flu (viral)?
 - Is it pneumonia (bacterial or viral)?
 - Is it both?

1. Bacterial Coinfections in Lung Tissue Specimens from Fatal Cases of 2009 Pandemic Influenza A (H1N1) — United States, May–August 2009: CDC MMWR, September 29, 2009; Vol. 58.
2. Metlay JP, et al. Diagnosis and Treatment of Adults with Community-acquired Pneumonia. An Official Clinical Practice Guideline of the American Thoracic Society and Infectious Diseases Society of America. Am J Respir Crit Care Med. 2019;200(7):e45-e67.

ATS/IDSA Clinical Practice Guidelines Recommend Influenza Testing in Patients with Pneumonia

In adults with CAP,
test for Influenza at the
time of diagnosis¹

 Strong
Recommendation

**Test for influenza with a
rapid influenza molecular assay,**
i.e., influenza nucleic acid amplification test;
preferrable to a rapid antigen test¹

In outpatients (including
emergency department)
and hospitalized patients²

**Clinicians should test for influenza in
patients with pneumonia²**

1. Metlay JP, et al. Diagnosis and Treatment of Adults with Community-acquired Pneumonia. An Official Clinical Practice Guideline of the American Thoracic Society and Infectious Diseases Society of America. Am J Respir Crit Care Med. 2019;200(7):e45-e67.
2. Uyeki TM, Bernstein HH, Bradley JS, et al. Clinical Practice Guidelines by the Infectious Diseases Society of America: 2018 Update on Diagnosis, Treatment, Chemoprophylaxis, and Institutional Outbreak Management of Seasonal Influenza. Clin Infect Dis. 2019 Mar 5;68(6):e1-e47.

CLINICAL, ECONOMIC AND OPERATIONAL

Improved Outcomes Associated With UAT

Impact of *Streptococcus pneumoniae* Urinary Antigen Testing in Patients With Community-Acquired Pneumonia Admitted Within a Large Academic Health System

Adam Greenfield,¹ Cassandra Marsh,¹ Justin Siegfried,¹ Ioannis Zacharioudakis,² Nabeela Ahmed,^{2,3} Arnold Decano,^{2,3} Maria E. Agüero-Rosenfeld,⁴ Kenneth Inglis,⁴ John Papadopoulos,¹ and Yanina Dubrovskaya^{1,2}

¹Department of Pharmacy, New York University Langone Health, New York, New York, USA, ²Division of Infectious Diseases, Department of Medicine, New York University Langone Health, New York, New York, USA, ³Department of Pharmacy, New York University Langone Hospital–Brooklyn, Brooklyn, New York, USA, and ⁴Department of Pathology, Grossman School of Medicine, New York University, New York, New York, USA

“We observed earlier de-escalation in the PUAT-positive group. ... due to discontinuation of atypical rather than anti-MRSA or antipseudomonal coverage.”

Shorter Duration of Antibiotic Therapy with Positive UAT Result

Table 3. Antimicrobial Days of Therapy During Entire Admission

Antimicrobial	Positive PUAT (n = 121)	Negative PUAT (n = 789)	PValue
Azithromycin	2 (1–3)	3 (1–4)	.024
Doxycycline	2 (1–3)	3 (2–4)	.027
Vancomycin	3 (1–4)	2 (2–4)	.908
Piperacillin-tazobactam	3 (2–6)	4 (3–7)	.053
Cefepime	1 (1–4)	1 (1–4)	.370
Ceftriaxone	4 (3–7)	2 (3–4)	.0005
Fluoroquinolone	2 (1–9)	2 (1–4)	.649
Linezolid	1 (1–2)	2 (1–8)	.272
Meropenem	3 (1–12)	5 (3–8)	.397
Ampicillin-sulbactam	1 (1–1)	1 (1–2)	.564
Broad-spectrum days of therapy			
Atypical coverage	2 (1–3) n = 103	3 (2–4) n = 722	.007
MRSA coverage	2 (1–4) n = 64	2 (2–4) n = 368	.625
<i>Pseudomonas aeruginosa</i> coverage	3 (2–5) n = 61	4 (2–6) n = 368	.315

Data are presented as median (interquartile range) unless otherwise stated; Antimicrobial exposure was determined throughout entire admission and patients may have received multiple agents.

Abbreviations: MRSA, methicillin-resistant *Staphylococcus aureus*; PUAT, pneumococcal urinary antigen test.

More Rapid Antibiotic De-escalation and Discontinuation with Positive UAT Result

Table 4. Comparison of De-escalation Between Pneumococcal Urinary Antigen Test–Positive and –Negative Groups

Characteristic	All Patients (N = 910)		PValue
	Positive PUAT (n = 121)	Negative PUAT (n = 789)	
Overall initial de-escalation	97/117 (82.9)	629/775 (81.2)	.746
Time to de-escalation from PUAT, d, median (IQR)	1 (0–2)	1 (1–2)	.01
Atypical coverage	n = 103	n = 722	
Discontinuation	80/103 (77.7)	509/722 (70.5)	.165
Within 24 h of PUAT	49/80 (61.3)	240/509 (47.2)	.026
Time to discontinuation, median (IQR)	1 (1–2)	2 (1–2)	.04
MRSA coverage	n = 64	n = 368	
Discontinuation	45/64 (70.3)	265/368 (72)	.898
Within 24 h of PUAT	24/45 (53.3)	127/265 (47.9)	.610
Time to discontinuation, d, median (IQR)	1 (1–2)	2 (1–2)	.131
<i>Pseudomonas aeruginosa</i> coverage	n = 61	n = 368	
De-escalation ^a	35/61 (57.4)	177/368 (48.1)	.228
Within 24 h of PUAT	20/35 (57.1)	99/177 (55.9)	.895
Time to de-escalation, d, median (IQR)	1 (1–2)	1 (1–2)	.621

Data are presented as No. (%) unless otherwise stated.

Abbreviations: IQR, interquartile range; MRSA, methicillin-resistant *Staphylococcus aureus*; PUAT, pneumococcal urinary antigen test.

^aDe-escalation defined as ≤ 3 days of therapy (discontinued within 3 days from initiation of antibiotic).

Improved Antibiotic Use and Length of Stay Achieved

Effect of a 3-Step Critical Pathway to Reduce Duration of Intravenous Antibiotic Therapy and Length of Stay in Community-Acquired Pneumonia A Randomized Controlled Trial

Jordi Carratalà, MD; Carolina Garcia-Vidal, MD; Lucía Ortega, MD; [et al](#)

OBJECTIVE CRITERIA USED FOR SWITCHING TO ORAL ANTIBIOTIC THERAPY

Pneumonia order set included:

- *Legionella* UAT
- *S. pneumoniae* UAT

RANDOMIZED TRIAL RESULTS:

- Median duration of IV antibiotic therapy **reduced 50%** (2.0 vs 4.0 days)
- Fewer patients experienced adverse drug reactions
- Median LOS **decreased 22%** (3.9 days vs 6.0 days)

LOS, Length of stay

Carratalà J, Garcia-Vidal C, Ortega L, et al. Effect of a 3-step critical pathway to reduce duration of intravenous antibiotic therapy and length of stay in community-acquired pneumonia: a randomized controlled trial. Arch Intern Med. 2012 Jun 25;172(12):922-8.

Key Takeaways for UAT


Guidelines Supporting UAT

ATS/IDSA CAP GUIDELINES¹

UAT recommended in:

Adults with **severe CAP** → *S. pneumoniae* and *Legionella*

Where indicated by **epidemiological factors** → *Legionella*



Evidence Supporting UAT | Reported Outcomes

Reduced **mortality**^{1,2,3}

- **57%** reduced odds of in-hospital mortality²
- **25%** reduced odds of 30-day mortality in UAT patients³

Improved
antibiotic stewardship
with positive UAT

- **125.9%** increase in narrowed therapy⁴
- **57%** of patients de-escalated antibiotic therapy¹
- **15%** of patients changed antibiotic therapy²

Shorter **Hospital
Length of Stay (LOS)**

- **57.8%** reduction¹

Fewer **adverse drug reactions**

Avoidance of ***Legionella* outbreaks**⁵

Why UAT for *S. pneumoniae* and *Legionella*?

CLINICAL RATIONALE

Guideline-concordant testing¹
Guide for antibiotic de-escalation^{2,3}
Demonstrated clinical utility and
associated with mortality outcomes²⁻⁶

UAT IS SUPPORTED BY THE ATS/IDSA,
ENHANCES ANTIBIOTIC STEWARDSHIP
AND IS ASSOCIATED WITH
IMPROVED CLINICAL OUTCOMES

LABORATORY RATIONALE

Non-invasive, ease of urine sample collection²
No instrument required, easy to use^{2,3}
Rapid results^{2,3}
Low cost per test^{2,3}

UAT IS A COST-EFFECTIVE, REQUIRES
NO INSTRUMENT TO MINIMIZE MOLECULAR
WORKFLOW DISRUPTION AND OPTIMIZE
UTILIZATION OF LIMITED RESOURCES

1. Metlay JP, et al. Am J Respir Crit Care Med. 2019;200(7):e45-e67.
2. West, et al. ERJ Open Res 2016; 2: 00011-2016.
3. Schimmel JJ, et al. Clin Infect Dis. 2020 Sep 12;71(6):1427-1434.
4. Costantini E, et al. Intern Emerg Med. 2016;11(7):929-940.
5. Uematsu H, et al. 2014;26(1):100-107.
6. Puri S, et al. 2020;17(2):533. Published 2020 Jan 15.

Summary

- High mortality rates, particularly *Legionella*
- *S. pneumoniae* is most likely bacterial cause of CAP
- Guidelines support the use of UAT
 - In severe CAP (*S. pneumoniae* and *Legionella*)
 - In cases of outbreak or recent travel (*Legionella*)
- Demonstrated clinical outcome improvements include
 - Reduction in in-hospital and 30-day mortality
 - De-escalation of broad-spectrum antibiotic treatment
- Shorter time to optimizing therapy associated with additional improved outcomes
 - Fewer adverse drug reactions
 - Shorter hospital length of stay
 - Reduced risk of *C. difficile* infection
 - Reduced cost of care

Thank You!

Appendix – Additional Citations

EVIDENCE SUPPORTING UAT – REPORTED OUTCOMES

1. West DM, et al. Pneumococcal urinary antigen test use in diagnosis and treatment of pneumonia in seven Utah hospitals. *ERJ Open Res.* 2016 Oct 19;2(4):00011-2016.
2. Costantini E, Allara E, Patrucco F, Faggiano F, Hamid F, Balbo PE. Adherence to guidelines for hospitalized community-acquired pneumonia over time and its impact on health outcomes and mortality. *Intern Emerg Med.* 2016;11(7):929-940.
3. Uematsu H, Hashimoto H, Iwamoto T, Horiguchi H, Yasunaga H. Impact of guideline-concordant microbiological testing on outcomes of pneumonia. *Int J Qual Health Care.* 2014;26(1):100-107.
4. Schimmel JJ, Haessler S, Imrey P, Lindenauer PK, Richter SS, Yu PC, Rothberg MB. Pneumococcal Urinary Antigen Testing in United States Hospitals: A Missed Opportunity for Antimicrobial Stewardship. *Clin Infect Dis.* 2020 Sep 12;71(6):1427-1434.
5. Puri S, Boudreaux-Kelly M, Walker JD, Clancy CJ, Decker BK. Clinical Presentation of Community-Acquired Legionella Pneumonia Identified by Universal Testing in an Endemic Area. *Int J Environ Res Public Health.* 2020;17(2):533.

Questions



Available CE Credit

P.A.C.E.[®]

Florida laboratory CE

Certificate of Attendance

To Receive Certificate of Attendance

After today's webinar:

- A certificate of attendance available for all attendees
- Evaluation form will appear automatically
- Must complete Eval to receive Certificate link via email
- For groups: Those logged in will receive Email from messenger@webex.com with link to evaluation. Forward email to colleagues who attended with you!!!
- Double-check email address

Joined Using a Mobile Device?

- Evaluation won't appear automatically, but...
- Watch for email with link to evaluation!



Recording

Within a few days following today's event, visit
<https://www.whitehatcom.com/abbott>



Know the Pneumonia: Urinary Antigen Testing in Patient Care

Live Event: Thursday, November 9, 2023 | 1:00 - 2:00 PM Eastern Time

P.A.C.E.® credit available until November 9, 2024

Florida Laboratory CE Credit available

RECORDING

SLIDES

Join this session for insights on improving the quality of care of patients with pneumonia. Experts will discuss the impact of this deadly infection and current gaps in identifying *Legionella* and *Streptococcus (S.) pneumoniae*. Practical applications of integrating urinary antigen testing (UAT) in patient scenarios will be explored, including an evidence-based review of mortality and antibiotic stewardship improvements achieved with UAT.

Presenter:



Janet E. Stout, PhD

Executive V.P. and Founder
Special Pathogens Laboratory
Pittsburgh, PA

Moderator/Speaker:



Norman Moore, PhD

Director, Infectious Diseases,
Medical Affairs, Rapid Diagnostics
Abbott
Scarborough, ME

The webinar will:

- Evaluate the impact of pneumonia and risks for severe illness
- Examine gaps and potential delays in the diagnosis of pneumonia and quality measures to enhance patient care
- Review the identification of *Legionella* and *S. pneumoniae* with UAT
- Discuss practical examples where UAT testing supported improved patient outcomes through timely identification

Know the Pneumonia: Urinary Antigen Testing in Patient Care

NOTE: If you have just viewed the archived recording of this webinar, you can access the evaluation using the link in the email you received after submitting the recording request form. Alternatively, you can access the evaluation for **12 months** after the live event at:

https://www.whitehatcom.com/Abbott_Evals/Pneumonia_110923/Testing_110923_eval.html