

# Deep Venous Thrombosis Pulmonary Embolism, D-dimer and Point-of-Care


Howard J. Kirchick, Ph.D.  
Director of Scientific Affairs  
Alere

# Objectives

- Review of Pathophysiology of DVT and PE
- Diagnosis
  - History & Physical examination
  - Imaging
  - Lab work
- D-dimer Tests
  - Latex
  - Immunometric
  - Specificity
  - Point of Care

# Venous Thromboembolism


- A blood clot, or thrombosis, develops abnormally in the blood vessel; usually the extremities.
- A deep vein thrombosis (DVT) forms primarily in the deep calf or thigh veins behind a valve.
  - May cause swelling if it persists
  - Most are relatively minor and go unnoticed
  - Pain occurs once extended along the vein and enters into thigh vein
- If DVT is not treated immediately, the blood clot may reach the lungs and cause a potentially fatal pulmonary embolism
- 90% of blood clots resulting in a PE stem from a DVT



# Venous Thromboembolism (VTE)


- 3rd most common cardiovascular disease
- Encompasses deep vein thrombosis (DVT) and pulmonary embolism (PE)

Ileo-femoral DVT



# Pulmonary Embolism (PE)

- Clots that travel through the venous system to reach and block a pulmonary vessel.
- If a clot reaches the pulmonary arteries, blood circulation is disturbed and subsequently gas exchange is partly hindered



# Partial List of Risk Factors


- Age >40 yr
- History of VTE
- Surgery/Trauma
- Prolonged immobilization
- Congestive heart failure
- Fracture of pelvis, femur or tibia
- Cancer
- Obesity
- Pregnancy or recent delivery
- Oral contraceptives/Estrogen therapy
- Inflammatory bowel disease
- Burns
- Genetic or acquired thrombophilia



# Clinical Symptoms of PE and/or DVT

- Shortness of breath 73%
- Chest Pain 66%
- Leg Pain or Swelling 33% (due to DVT)
- Cough 43%, sometimes with blood 15%
- Tachycardia
- Dizziness
- Syncope
- Tachypnea
- Crackles
- Jugular venous distention
- Fever
- EKG changes

**#1 method for diagnosis:**






# Issues in Diagnosing Patients with SOB

- Differential diagnoses
  - PE
  - Myocardial Infarction
  - Congestive Heart Failure
  - Pneumonia
  - COPD
  - Cardiac Tamponade
- Diagnostic testing
  - Cardiac markers
  - D- dimer
  - CBC, chemistry, lipid panel
  - EKG
  - CXR
  - VQ scan/CT scan
  - Cultures
  - Echocardiogram
  - Stress test
  - Left/Right Heart Catheterization
  - Pulmonary Function Test

# Other Causes of D-dimer Elevation

- DIC
  - AMI
  - Atherosclerosis
  - Trauma
  - Hepatic disease
  - Sepsis
  - Surgery
  - Infection
  - Pregnancy
  - Inflammation
  - Age
  - Cancer
  - Thrombolytic therapy
- Hence a positive test does not prove the existence of DVT/PE**

# Current practice in PE diagnosis?



## General disadvantages:

1. Instrument and skilled staff have to be available
2. Potential of renal damage as a result of imaging dye administration

# So What's The Problem?

- The clinical presentation of both DVT and PE may be misinterpreted, subtle or asymptomatic
- Radiologic studies are expensive, subjective and often non-diagnostic, potentially harmful to the patient and not always readily available
- Need a simple, fast, inexpensive test that is highly sensitive and preferably specific

# Challenges Associated With D-dimer

- Not specific to a disease, detects breakdown of clot
  - May encounter false positives for PE/DVT
  - Value is in ability to reduce further evaluation of patients with a negative D-dimer
- No standard for D-dimer; results vary, correlation is difficult
  - Latex agglutination subjective and has 80% sensitivity, vs. sandwich immunoassays (ELISA and FIA) assays with nearly 100% sensitivity

# How is D-dimer Being Used ?

## D-dimer often used inappropriately

- No pretest probability assessment
- Blanket test of all chest pain patients
  - Overuse diminishes the value of the test
    - Lowers specificity / increases FP rate, decreases clinician/lab confidence in test!
  - Shortens life of the scanner
    - With D-dimer screening, the positive rate of CT scans for PE is 11%\* - 15%\*\*
    - Without D-dimer screening, the positive rate is 5%\*\* - 8%\*
  - Irradiates patients
    - The radiation from one chest CT  $\geq$  40 chest x-rays
- D-dimer TATs insufficient to make rapid clinical decisions for imaging
  - Ordering D-dimer, but sending concurrently for imaging, if available

\* Kline, et al, *Annals of Emergency Medicine*, Nov 2004

\*\* Night Radiologist et al. Sharp Hospital. Unpublished

# Appropriate Use of D-dimer

## American College of Emergency Physicians Clinical Policy Statement

- In most cases, low probability patients are candidates
- Screen patients with a Pre-test Probability Score (Wells, Hamilton, Charlotte, Geneva, etc)
- Use in out-patient population
  - Hospitalized, pregnant, post surgical patients will likely be elevated due to other clinical conditions/risk factors
- When used appropriately, D-dimer assists in reducing the number of patients requiring CT scans
- Physician education will be VERY useful
  - Use on low probability patients that would otherwise be sent for imaging/scanning as part of a PE or DVT workup

# ACEP Clinical Policy

- In patients with low pre-test probability the following can be used to exclude PE:
  - Negative quantitative D-dimer
  - Negative whole blood qualitative D-dimer AND Wells' score < 2

**Low Probability** <2.0  
*(3.6% Risk)*

**Moderate** 2.0 – 6.0  
*(20.5% Risk)*

**High Probability** >6.0  
*(66.7% Risk)*


Clinical Characteristics	Score
Clinical signs and symptoms of DVT	3
PE likely or more likely than alternative diagnosis	3
Heart rate greater than 100 beats/min	1.5
Immobilization (bedrest ≥ 30 days) or surgery in the previous 4 weeks	1.5
Previous DVT/PE	1.5
Hemoptysis	1.0
Malignancy (Receiving treatment, treated in the last 6 months, or palliative care)	1.0

# Strategy For Diagnosis of PE



Note: Major role for D-dimer is the low risk ED or outpatient

Note: High risk or inpatient: little role for D-dimer: But remember DIC





# **Implementation of a Rapid Whole Blood D-Dimer Test in the Emergency Department**

Lewandrowski *et al.*, *Am J Clin Pathol* 2009;132:326-331



# Rapid Whole Blood Test in the ED

- Study Objectives – To Assess:
  - ED length of stay pre- and post-implementation of POC D-dimer;
  - Admission and discharge rates pre- and post-implementation of POC D-dimer; and
  - Utilization of imaging test rates pre- and post-implementation of POC D-dimer.

# Methods and Materials

- 252 patients pre-implementation and 211 patients post-implementation were evaluated for:
  - test results, turnaround times, and test volumes
  - ED LOS
  - patient chart reviews

# Results

- Following implementation of the rapid D-dimer test the total test turnaround time (from blood draw to availability of the test result) decreased from approximately 2 hours (central laboratory, depending on the shift and time of day) to 25 minutes, representing an approximately 79% decrease.

# Results

- The volume of D-dimer tests requested by the ED increased from a mean of 127 per month before implementation of the rapid D-dimer test to a mean of 154 tests per month (a 21.3% increase;  $P = 0.037$ ), reflecting increased utilization.
- Some of this increase can be explained by an approximately 6% increase in ED visits during the study period (daily average of 221 before to 235 after implementation).

# Results

## Rates of Hospital Admission, Discharge, and Admit to Observe for Patients Before and After Implementation of the Rapid Whole Blood D-Dimer Test in the Emergency Department

	Before Implementation	After Implementation
Admitted (%)	36.5	22.7
Discharged (%)	42.9	50.2
Admit to observe (%)	20.6	27.0

- The difference pre- and post-implementation was significant ( $P = 0.005$ ), indicating that the availability of the rapid test may have influenced patient disposition decisions.

# Results

## Rates of Follow-up Radiologic Testing\* Before and After Implementation of the Rapid Whole Blood D-Dimer Test in the Emergency Department

Radiologic Study	Before Implementation	After Implementation
No (%)	60.3	61.1
Yes (%)	39.7	38.9

\* Venous ultrasound, lung scan, or computed tomography.

- There was no statistical difference in imaging rates, i.e., implementation of POC D-dimer did not increase imaging.

# Key Points

- The POC test performed as well as the Lab test while producing
  - A significantly shorter ED LOS
  - Fewer admissions
  - No change in the rate of imaging
- These changes should result in decreased costs.



# **Diagnostic Accuracy and User-Friendliness of 5 Point-of-Care D-Dimer Tests for the exclusion of Deep Vein Thrombosis**

Geersing *et al. Clin Chem* 2010;56:1758-66

.



# DVT Rule Out Study

- Study Objectives – To Assess:
  - Accuracy of five POC D-dimer tests in the assessment of deep vein thrombosis (DVT);
  - Nursing opinion on the ease of use of the tests.

# Methods

- Protocol
  - 577 patients suspected of having a DVT were evaluated using:
    - Vidas (bioMerieux)
    - Pathfast (Mitsubishi)
    - Triage (Alere)
    - Cardiac (Roche)
    - Clearview Simplify (Alere)
- Ease of Use
  - Twenty nurses completed questionnaires dealing the ease of use of the five test methods.

# Results

- Accuracy:
  - Differences in the calculated sensitivities and specificities and negative predictive values were largely a result of differences in the cutoffs used.
  - There were no significant differences in the area under the curves (AUC) for the ROC analyses.

**Table 2. Diagnostic accuracy measures of the 5 point-of-care D-dimer tests for the exclusion of DVT (N = 577).**

D-dimer assay	Cutoff value	Sensitivity (95% CI)	Specificity (95% CI)	LR <sup>-</sup> (95% CI)	NPV (95% CI)	AUC (95% CI)
Vidas	500 µg/L FEU	0.99 (0.96–1.0)	0.42 (0.37–0.46)	0.03 (0.01–0.24)	99.5 (98–100)	0.89 (0.82–0.97)
Pathfast	570 µg/L FEU	0.98 (0.94–1.0)	0.39 (0.35–0.44)	0.06 (0.02–0.29)	99.2 (98–100)	0.89 (0.81–0.97)
Triage	196 µg/L FEU	0.97 (0.93–1.0)	0.48 (0.44–0.53)	0.06 (0.03–0.23)	99.2 (98–100)	0.88 (0.80–0.97)
Cardiac	500 µg/L FEU	0.94 (0.88–0.99)	0.62 (0.58–0.67)	0.10 (0.07–0.25)	98.6 (97–100)	0.88 (0.80–0.96)
Clearview Simplify	80 µg/L	0.91 (0.85–0.98)	0.64 (0.60–0.69)	0.14 (0.09–0.29)	98.1 (97–100)	NA

<sup>a</sup> LR<sup>-</sup>, likelihood ratio negative; NPV, negative predictive value; NA, not applicable because it is a qualitative assay.


# Results

- Ease of use:
  - “On the basis of user-friendliness the Cardiac and Triage devices may be preferred for use in an emergency department setting or small primary care clinics. In addition, these tests produce a D-dimer test result within 15 min and can provide measurements of other (cardiovascular) biomarkers ...”
  - “The Clearview Simplify test is also user-friendly for primary care, because it is easily portable, requires no analyzer, can be performed on capillary whole blood, and requires no calibration.”
  - The drawback to the Clearview Simplify is that the interpretation of test results was rated more difficult because of the subjective nature of the reading.

# D-dimer tests-Choices and Challenges

- Latex Agglutination
  - Qualitative or semi-quantitative
  - Relatively insensitive
- Whole Blood Agglutination
  - More sensitive than latex
  - Somewhat subjective
- Turbidimetric
  - Much less subjective
- ELISA/Immunoassay
  - High analytical sensitivity
  - Historically laborious
- Standardization challenges
  - No recognized standard
  - Purified D-dimer from plasmin digested clots (ng/mL)
  - D-dimer from totally lysed fibrin clots (Fibrin Equivalent Units, FEU)
- Patient variables

# Tests of Fibrinolysis




- Fibrinogen
- Platelet count
- Fibrin degradation products
  - FpA
  - FpB
  - Fragment D
  - Fragment E
- D-dimer (“cross-linked” fibrin degradation product)

■ Only D-dimer is useful for DVT and PE

# How Were D-Dimers Measured?

- Latex Agglutination:
  - Big clumps that are visible to the naked eye



- Turbidimetric assays:
  - Big clumps that scatter light – the less light detected, the more analyte is present

# Turbidimetric Assays

- Shine a light on one side and measure the light coming through on the other side




# ELISA / EIA

- Enzyme-Linked Immunosorbent Assay
  - Synonymous with Enzyme Immunoassay (EIA)
  - 1st ELISAs were run in microtiter plates (aka ELISA plates)
- Member of a class of immunoassays (Immunometric or “Sandwich”)
  - All involve capturing the analyte
  - All involve measuring captured analyte using a form of signal generator
    - EIA uses an enzyme-labeled antibody to convert an “invisible” molecule into a “visible” molecule
    - FIA (Fluorescence immunoassays, or IFA, immunofluorometric assay) are similar to EIA except that they use a fluorescent-labeled antibody as the signal
    - ELFA-Enzyme linked fluorescent immunoassay
  - FIA can be just as sensitive as EIA (e.g., TnI or BNP)



# “Sandwich” Immunoassay

- Typically used for analytes with multiple epitopes (Cardiac Markers, D-dimer, Microbiology).





# Assays Compared

- Shown below is a POC immunometric (sandwich) assay versus an immunometric Lab assay (left panel) and the POC assay versus a turbidimetric Lab assay (right panel).



# Value of D-dimer Antibody Specificity


- False positives reduce the value of D-dimer and increase clinician and lab frustration.
- Tests with high affinity antibodies for D-dimer reduce false positives.



- The 3B6 monoclonal antibody offers high specificity due to its affinity to the cross-linking epitope (recognition site) of D-dimer.

# Review of 78 DVT/PE Studies


- 78 prospective clinical studies investigated the use of D-dimer for the exclusion of acute VTE and PE
- The specificity the 3B6-based whole blood assay was identified as clinically and statistically superior to the rapid ELISA and automated latex immunoassay methods for acute DVT and PE.



# Fibrin Assay Comparison Trial (FACT)

- Study Findings:
  - The main reason for differences between D-dimer assays is due to differences in antibody specificity
  - Assays displaying cross-reactivity with non-cross linked fibrinogen and fibrin derivatives will show falsely high
    - Diagnostica Stago assays (MAbs 8D2, 2.1.16) showed greater than 30% cross-reactivity
  - Assays using 3B6 antibodies were identified as the most specific for D-dimer. 3B6 assays had the least false positives.

## Distinguish from other FDPs



- FpA
- FpB
- D-dimer
- Fragment D
- Fragment E

Plasmin-derived FDPs may be detected in addition to D-dimer, resulting in an erroneously elevated result

# Distinguish from other FDPs



- False positives reduce the value of D-dimer and increase clinician and lab frustration
- Tests with high affinity antibodies for D-dimer reduce false positives
- The 3B6 monoclonal antibody offers high specificity due to its affinity to the cross-linking epitope of D-dimer

# Is D-dimer useful for DIC?

Diagnostic algorithm for the diagnosis of overt disseminated intravascular coagulation.

- Risk assessment: Does the patient have an underlying disorder known to be associated with overt DIC?
- If yes, proceed. If no, do not use this algorithm;
  1. Order global coagulation tests (platelet count, PT, fibrinogen, soluble fibrin monomers, or fibrin degradation products).

# Is D-dimer useful for DIC?

## 2. Score global coagulation test results:

- Platelet count
  - ( $>100 \times 10^9/L = 0$ ,  $<100 \times 10^9/L = 1$ ,  $<50 \times 10^9/L = 2$ )
- Elevated fibrin-related marker (e.g. soluble fibrin monomers/fibrin degradation products - **D-dimer**)
  - Historical abnormal D-dimers can be split into tertiles
    - no increase scores a 0; mild increase (1<sup>st</sup> tertile) scores a 1; moderate increase (2<sup>nd</sup> tertile) scores a 2; strong increase (3<sup>rd</sup> tertile) scores a 3.
- Prolonged prothrombin time
  - ( $<3 \text{ s} = 0$ ,  $>3 \text{ but } <6 \text{ s} = 1$ ,  $>6 \text{ s} = 2$ )
- Fibrinogen level
  - ( $>1.0 \text{ g/l} = 0$ ,  $<1.0 \text{ g/l} = 1$ )

# Is D-dimer useful for DIC?

## 3. Calculate score.

- If  $\geq 5$ : compatible with overt DIC; repeat scoring daily.
- If  $< 5$ : suggestive (not affirmative) for non-overt DIC; repeat next day.

# D-dimer, conclusions

- Most appropriate for ED patients as hospitalized will usually have elevated levels
- When used appropriately, D-dimer is a useful tool for ruling out venous thromboembolism and reducing costs and adverse outcomes that result from unnecessary imaging studies.
- D-dimer can be used on in-patients to help assess disseminated intravascular coagulation (DIC).

# Questions????

